


Potentiometric Titrations calculations.



The original E_{obs} is the value recorded with the silver wire and the reference electrode in the 0.070 M $AgNO_3$ solution and is used to calculate E_{ref} .

*Note there are three E_{obs} values to be used in the calculations at various times.

1. Use your titration curve to determine the two equivalence points and to find the I^- and Cl^- concentrations in your unknown samples. Note: The E_{obs} values to be used in K_{sp} calculations are also determined at these points. The equivalence points can be more accurately determined by the first and/or second derivatives.

2. $E^o = 0.80$ V. E_{ind} must be determined using the formula: $E_{ind} = E^o - \frac{RT}{nF} \ln[Q]$ Where $Q = 1/[Ag^+]$ ($[Ag^+] = 0.070$ M in your case), $R = 8.314$ J/mol·K, $n = 1$, $F = 9.65 \times 10^4$ C/mole of e^- .

3. Use E_{ind} and your E_{obs} value (the reading you got with the silver wire and the reference electrode and the standard 0.070 M $AgNO_3$ solution) in the formula: $E_{ind} - E_{obs} = E_{ref}$ to calculate E_{ref} .

4. Now use the formula (Nerst Equation): $[Ag^+] = e^{\frac{nF}{RT}(E_{obs} + E_{ref} - E^o)}$ to determine $[Ag^+]$. The E_{obs} in this equation is the potential at the equivalence point and E_{ref} is the value determined in number 3.

5. Then, since $[Ag^+] = [I^-]$ for AgI at the equivalence point, the K_{sp} of AgI is equal to $[Ag^+]^2$.

6. Repeat the calculations done in 4 and 5 for $AgCl$ using the E_{obs} from the 2nd equivalence point to determine the experimental K_{sp} of $AgCl$.

7. Look up the accepted values of the two K_{sp} 's and compare them to you values.

8. When finished, you should be reporting four items: $[Cl^{1-}]$, $[I^{1-}]$, K_{sp} of $AgCl$, and K_{sp} of AgI .

